

Novel SPE Extraction Method for Sensitive and High Throughput Quantitative Analysis of Phosphorothioate Oligonucleotides in Human Plasma Using LC-MS/MS

Yuhuan Ji; Laixin Wang; Qian Liu; Qian Li; Ji Chengjie
NovaBioAssays, LLC, Woburn, MA

Overview

A high throughput SPE extraction method using a mechanism that is different from conventional reverse phase and ion exchange mechanism has been developed.


Introduction

Oligonucleotides and RNA are gaining renewed confidence as a new class of drugs. Many bioanalytical methods have been developed to support drug research and development of oligonucleotide therapeutics. LC-MS/MS is one of the most popular techniques due to its unprecedented specificity; however, the sensitivity for most LC-MS/MS assays is still not as good as qPCR and ligand binding assays for some large oligonucleotides. To improve sensitivity and throughput of the LC-MS/MS assays, we have successfully developed a novel SPE extraction method which produces high quality extracts from human plasma samples and requires no drying down step.

Methods

A 20-mer phosphorothioate DNA oligonucleotide (NBA1520, Mw = 6387.1 Daltons) and a 14-mer analog oligonucleotide (NBA1514, Mw = 4437.6 Daltons) were used as a model analyte and internal standard, respectively. To optimize the SPE conditions, the loading buffer, washing buffer, and elution buffer were screened and optimized by extracting 200 µL of 5,000 ng/mL NBA1520 in human plasma using the Biotage EVOLUTE EXPRESS WAX 30 mg Fixed Well plate (Part No.: 604-0030-PX01). The eluents were analyzed directly on a Thermo DNAPac C18 column (2x50 mm, 4 µm) with a gradient of HFIP/TEA buffered water and methanol. The standard curve (0.5 ng/mL-500 ng/ml) and QC samples were quantitatively analyzed by SCIEX API5000 (MRM-MS) and Thermo Q Exactive plus (targeted-SIM) under negative ion electrospray mode.

Experimental Workflow



Conclusions

- ❖ A high throughput SPE extraction method have successfully developed to purify phosphorothioate DNA oligonucleotide from human plasma
- ❖ LC-HRAM on Q-Exactive plus can achieve comparable sensitivity as LC-MS/MS on triple quads for same analytes

Results

